




## DEMS System

► Differential Electrochemical Mass Spectrometry



# Real Time Gas Analysis for Advanced Energy and Electrocatalytic Research – DEMS system

Differential Electrochemical Mass Spectrometry is a technique that allows unique online analysis of gaseous and volatile products from electrochemical reactions. The Hiden DEMS systems combines advanced electrochemical half cells with excellent performance in mass spectrometry to achieve outstanding analytical capabilities.

Designed for applications such as fuel cell diagnostics, CO<sub>2</sub> reduction, catalyst evaluation, corrosion studies and battery research. It provides scientists with accurate, immediate data on reaction pathways and product selectivity.



## Key Features

- ▶ Versatile Mass Spectrometer system compatible with multiple cells and inlet types
- ▶ Near-instantaneous time response – direct electrode-to-MS connection delivers rapid, quantitative product detection
- ▶ Mass scanning and time/intensity trend monitoring of multiple species
- ▶ Deposition compatible working electrodes
- ▶ Data synchronisation within a single intuitive software package
- ▶ A turnkey solution – everything you need to start your experiments immediately

## System Configuration & Options

ITEM	DESCRIPTION	PARTCODE
SYSTEM	HPR-40 DEMS bench-top gas analysis system for electrochemistry, including ECL-Static and ECL-Insight cells and Hiden HAL 201 RC mass spectrometer with Faraday/Electron Multiplier detector. Mass range 200 amu.	305250
OPTIONS & ACCESSORIES	Extended mass range. 300 amu mass range (in place of standard 200 amu mass range).	305021
	<b>ECL-Probe</b> – Extendable probe, allowing dissolved gas to be measured in custom cells.	303443
	<b>Potentiostat Integration</b> – Cables available for connection to most Potentiostats, allowing trigger start and real time integration of Potential and Current into Hiden software.	270220/270221
GAS INLET OPTIONS	<b>QIC inlet</b> - heated capillary inlet for sampling gases and vapours at atmospheric pressure.	303562
	Microflow inlet, flow rate from 12 µl/min, unheated.	303452
SPARES KIT	<b>Recommended spares kit</b>	
	▶ Twin filament, oxide coated iridium	201200
	▶ Filament kit	201600
	▶ ECL-Static	3035515
	▶ ECL-Insight	3035517

# ECL-Series





## ECL-Static

ECL-Static cell is optimised for ease of use and flexible experimentation setup. Its four ports accommodate counter and reference electrodes as well as allowing for electrolyte flow, should a flow configuration be required. Depending on the application, the vitreous carbon working electrode can be used as an inert electrode or an electrocatalyst can be applied to the active surface using various deposition techniques. The cell is supplied 'Research Ready' while still offering flexibility to customise electrodes and configurations.

## ECL-Insight

ECL-Insight cell is designed for high-performance, real-time electrochemical mass spectrometry, offering direct electrode-to-MS connection for rapid, quantitative detection. With a gold sputter-coated membrane, dual-chamber non-metallic design, and bubble-free operation at currents up to  $100 \text{ mA/cm}^2$ , it combines exceptional collection efficiency, ultra-fast voltage scans up to  $100 \text{ mV/s}$ , and high sensitivity down to  $\sim 50 \text{ } \mu\text{A/cm}^2$ . This makes it the ideal choice for studies in  $\text{CO}_2$  reduction, hydrogen evolution, and isotope labelling.



## ECL-Probe

The ECL-Probe enables direct dissolved gas analysis in both standard and custom electrochemical cells. Its non-metallic, immersion-ready construction ensures compatibility without electrochemical interference, while a specialised hydrophobic support delivers an enriched gas signal. Ready-to-connect with the HPR-40 DEMS, the ECL-Probe offers seamless integration and reliable performance.

# Example data




Figure 1: Real-time hydrogen generation during water electrolysis using ECL-Insight, showing the clear synchronisation of the hydrogen signal with applied potential and current in QGA software.






Figure 2: CO<sub>2</sub> reduction activity of Cu measured during cyclic voltammetry at a scan rate of 1 mV/s in 0.1 M CsHCO<sub>3</sub> at a flow rate of 85 mL/min.

# Technical data


**ECL-Insight electrode view**




**ECL-Insight electrolyte view**



**ECL-Static electrode view**



**ECL-Static electrolyte view**



Mass ranges, amu:

1-200 / 1-300 amu

Sensitivity:

100% to 100 ppb subject to spectral interference

Speed:

Up to 650 measurements/second

Response time:

<100 ms

Software:

MASsoft Professional

Windows compatible

QGAssoft

Interface:

Ethernet/USB/Serial (RS-232) connections

Detector:

Dual Faraday/Channeltron Electron Multiplier

Analogue input:

8x (optional)/16 bit

Analogue output:

8x (optional)/14 bit

Digital input:

8x

Digital output:

8x, 24 V

Dimensions (L x W x H), mm:

495 x 538 x 382 mm

Weight, kg:

Typically 31 kg and external scroll pump 26 kg

Power requirement:

110/220/240 V AC, 50/60 Hz, 1.2 kVA

# Hidden**APPLICATIONS**

Hidden's quadrupole mass spectrometer systems address a broad application range in:

## **GAS ANALYSIS**

- ▶ Dynamic measurement of reaction gas streams
- ▶ Catalysis and thermal analysis
- ▶ Molecular beam studies
- ▶ Dissolved species probes
- ▶ Fermentation, environmental and ecological studies




## **SURFACE ANALYSIS**

- ▶ UHV TPD/TDS
- ▶ ToF qSIMS and SIMS analysers
- ▶ End point detection in ion beam etch
- ▶ Elemental imaging – 3D mapping
- ▶ SIMS system with simultaneous dual polarity analysis



## **PLASMA DIAGNOSTICS**

- ▶ Plasma source characterisation
- ▶ Etch and deposition process reaction kinetic studies
- ▶ Analysis of neutral and radical species



## **VACUUM ANALYSIS**

- ▶ Partial pressure measurement and control of process gases
- ▶ Reactive sputter process control
- ▶ Vacuum diagnostics
- ▶ Vacuum coating process monitoring



# **HIDEN**

**ANALYTICAL**

Hidden Analytical Ltd.  
420 Europa Boulevard  
Warrington WA5 7UN England

**T** +44 [0] 1925 445 225

**E** [info@hidden.co.uk](mailto:info@hidden.co.uk)

**W** [www.HiddenAnalytical.com](http://www.HiddenAnalytical.com)



### **Sales Offices:**

We have sales offices situated around the globe. **Visit our website for further information.**